



MGMR[®] Versatile Nano-carrier for Non-viral Cellular Transfection

Aaron Tasset, Kevin Castillo, Milos Marinkovic

Biopact, LLC

12310 Trail Driver St. 78737


Austin, Texas, USA



MGMR[®] Supports T-cell Proliferation


Dose-response Testing for Proliferation of Human T Lymphocyte-like cells (Jurkat)

- MGMR does not reduce cell proliferation relative to control
- (Growth Media, No Treatment)
- All groups indicated high viability
- Certain dosing-ranges of MGMR may improve cell proliferation and increase yield
- Protocol
 - Groups incubated with various concentrations of MGMR for 3 and 7 Days
 - Control = Standard Growth Media (RPMI 1640 + 15% FBS)
 - Cellular viability and quantification assessed by trypan exclusion assay
 - Stats: 95% CI



Endocytotic Trafficking



- MGMR is transported into the cell by active transport
 - Unlike other transport systems, MGMR does not damage the cell's external casing (plasma membrane) in order to enter inside
- MGMR is initially encapsulated into larger compartments called vesicles - that traffic into the interior of the cell
 - MGMR is cargo agnostic - it is not prevented from getting inside cells by the size of its molecular cargo
- Once internalized, MGMR is released from vesicles into the intracellular environment
 - MGMR can deliver cargo to many different locations inside the cell



MGMR[®] Peptide Delivery


MGMR demonstrated to enhance cell penetration of KLA—a cell “impermeable” 1.5kD protein^{1,2}

- Proteins are powerful cell-regulators but very challenging for intracellular delivery
 - Do not efficiently cross the cell's membrane
 - Protein activity is very sensitive and easily diminished
- In this case study, MGMR binds to KLA with high affinity, transports it across the plasma membrane, and releases it—triggering apoptosis
 - MGMR+KLA peptide produced absolute apoptosis
 - MGMR alone has no effect (high tolerability)
 - KLA alone has no effect (benign externally)




¹Zurita et al. Cancer Research 64, 435-439, 2004 ²Yang et al. Methods Mol. Biol. 1266, 29-53, 2015

MGMR[®] siRNA Delivery




More Dark Space = More Knockdown of Green Fluorescent Protein



- Unlike detergent-based delivery systems, MGMR is non-cytotoxic and compatible across many types of molecular cargo


MGMR[®] mRNA Delivery

mRNA Delivery to T Lymphocytes (Jurkat)





• 72-hour incubation • 1x = 500ng/mL • n = 3

mRNA Delivery to Monocyte-like Cells (U937)




• 72-hour incubation • 1x = 500ng/mL • n = 3

MGMR[®] Nuclear Staining in Mesenchymal Stem Cells



MGMR[®] Delivery DNA-binding Enzymes in T-cells




Cells Not Treated with DNA-binding Enzymes Do Not Fluoresce

MGMR Transports DNA-binding Enzymes into Cells

MGMR[®] Endonuclease-loading for Sequential DNA Cleavage

- Protocol**
- Scal+MGMR incubated with 1ug plasmid for 15 min at 37°C, followed by HincII+MGMR for additional 15 min at 37°C
 - Electrophoresis (70 min) - 1% agarose
- Results**
- Individual nucleases loaded on MGMR can cleave pDNA in series
 - MGMR preserves nuclease activity and demonstrates loading density-dependent pDNA cleavage
 - Sequential addition of “free” nucleases without MGMR produces very low-efficiency pDNA cleavage
 - Formulation 2 demonstrated greatest efficiency for serial plasmid cleavage relative to Formulation 1 and sequentially added free nucleases



MGMR[®] Next-generation Transfection Technology

| Cargo | MGMR™ | Viral Vectors | Polymers | Liposomes |
|---------------------|-------|---------------|----------|-----------|
| Small Molecules | + | - | + | 0 |
| Peptides | ++ | + | + | 0 |
| Nucleic Acids | ++ | ++ | + | 0 |
| Proteins | ++ | + | + | 0 |
| Lipopolysaccharides | ++ | - | + | 0 |
| Log(P) Independent | + | + | - | - |
| Properties | | | | |
| Non-immunogenic | + | - | + | + |
| Non-GMO | + | - | + | + |
| Simple Optimization | ++ | - | 0 | - |
| Scalable | ++ | + | + | - |
| Eukaryotes | + | - | 0 | 0 |
| Prokaryotes | + | - | 0 | - |
| Plants | + | - | 0 | - |

Key

- ++ Advantage
- + Compatible
- 0 Limited Compatibility
- Non-Compatible

How can MGMR[®] improve cell engineering?

- Eliminates** viral-mediated transfection
 - No risk of insertional mutagenesis
 - No viral toxicity and immunogenicity
 - Non-replicative
- Streamlines** manufacturing and quality processes
 - Non-biological transfection technology reduces variation
 - Readily extractable from final cell product
 - Eliminates need for new vector development – Single transfection vehicle for all cellular gene-editing systems
- Versatile** loading profile effective for intracellular delivery of:
 - Gene-editing complexes
 - Oligonucleotides
 - Peptides/Proteins
- Improves** cell yield
 - Reduced toxicity vs. viral or detergent-mediated transfection
 - Improved cell proliferation

MGMR[®] Discover its Value in Your Process

Seamlessly Combine MGMR[®] with your Technology in 3 Phases

- Send Biopact Your Payload Molecule**
 - Provide Biopact with your gene-editing technology:
 - CRISPR
 - TALEN
 - ZFN
 - Biopact loads MGMR with your payload at **No Cost**
- Biopact Evaluates Performance**
 - Biopact will determine:
 - Loading efficiency on MGMR
 - Transfection efficiency in your target cell
 - Biopact reports performance of MGMR to you at **No Cost**
- Put MGMR[®] to the Test**
 - Order MGMR loaded with your molecule (\$10k-\$20k)
 - Readily sufficient for *in vitro*, *ex vivo* or *in vivo* testing
 - Evaluate MGMR loaded with your payload in your own studies

biopact™